



Precice Bone Transport

The all-internal solution designed to treat segmental bone defects up to 10cm

Due to the limitations of current solutions in the market, the need for a reliable and simple approach for bone transport procedures is prevalent. The Precice® Bone Transport System delivers strength and stability for better patient satisfaction with this challenging procedure.¹

Less Invasive Solution


When compared to alternative solutions, the system provides a less invasive treatment for segmental bone defects.¹

All-internal device: Ability to achieve distraction osteogenesis through the intramedullary canal

Stability for a Vast Range of Defects

Intentionally designed to address a vast range of segmental defects, up to 10cm.

Built-in dual slot: Provides stability to support the transport of more complex segmental defects



Reduced Complexity for the Patient and Surgeon

When compared to alternative solutions, the system offers a simplistic approach to the treatment of complex segmental bone defects.

Simplistic solution: All-internal device provides a more simple solution for both the surgeon and the patient

Patient preferred: Patients prefer the Precice internal device with respect to overall satisfaction, reduced pain, ease of physical therapy and better cosmetic appearance¹

NuVasive Specialized Orthopedics, Inc.
101 Enterprise, Suite 100, Aliso Viejo, CA 92656 USA
+1 949.837.3600

Reference

1. Landge V, Shabtai L, Gesheff M, et al. Patient satisfaction after limb lengthening with internal and external devices. *J of Surg Orthop Advanc* 2015;24(3):174-9.

©2018. NuVasive, Inc. All rights reserved. 9501896 A

CE 0297

nuvative.com

System Offerings

Available Diameters

- 10.0mm
- 11.5mm
- 13.0mm

Available Styles

- Antegrade Femur Piriformis Straight
- Antegrade Femur Trochanter 10° Bend
- Retrograde Femur 10° Bend
- Retrograde Femur Straight
- Antegrade Tibia 10° Bend (10.0 and 11.5mm only)

Length	Stroke	Slot
280mm	60mm	Single: 60mm
300mm	70mm	Single: 70mm
320mm	80mm	Dual: 70 x 10mm
340mm	90mm	Dual: 70 x 20mm
360mm	100mm	Dual: 70 x 30mm
380mm	100mm	Dual: 70 x 30mm
400mm	100mm	Dual: 70 x 30mm

EC REP NuVasive UK Ltd.
Suite B, Ground Floor, Caspian House,
The Waterfront, Elstree, Herts WD6 3BS UK
+44 (0) 208.238.7850